

## Summary:

Modern medical devices include increasingly sophisticated and delicate components. This test is designed to evaluate the effect of an EOGas and Anprolene sterilization cycle on the functionality of a laptop computer. An abnormally high quantity of gas was used for each cycle to create a worst case. Back to back Anprolene and EOGas cycles were performed on the same laptop. NOTE: Battery was removed and wrapped separately prior to each cycle.

## Materials:



1. Dell Inspiron Laptop Model: PP05L
2. AN74i Sterilizer
3. EOGas 333 Sterilizer
4. AN71 Bag
5. 2 AN79 Ampoules
6. Biological Indicators
7. AN87/1087 Dosimeters
8. AN1005 Bag
9. 2 AN1005 Cartridges
10. CSR Wrap
11. Humidichip
12. Humiditube

## Method:

### Anprolene Cycle

1. The battery from the Dell computer was removed and wrapped in CSR wrap.
2. The laptop, wrapped battery, AN87 Dosimeter, Biological Indicator, Humidichip within a Humiditube and two (2) AN79 Ampoules were placed within an AN71 Bag. NOTE: The AN71 bag is normally used with a 5 gram ampoule of EO. Two AN79 ampoules represent 35 grams of EO, or 700% more gas than a normal cycle.
3. The bag was loaded into an AN74i Sterilizer and a 12 hour sterilization cycle was completed.
4. Upon completion of the sterilization cycle the biological indicator was incubated for seven (7) days and the computer was tested for functionality.

### EOGas Cycle

1. The battery from the Dell computer was removed and wrapped in CSR wrap.
2. The laptop, wrapped battery, AN1087 Dosimeter, Biological Indicator, Humidichip within a Humiditube and two (2) AN1005 Cartridges were placed within an AN1005 Bag. NOTE: The AN1005 bag is normally used with a single 1005 (4 gram) cartridge of EO. Two AN1005 cartridges represent 8 grams of EO, or 200% more gas than a normal cycle.
3. The bag was loaded into an EOGas 333 Sterilizer and a 16 hour sterilization cycle was completed.
4. Upon completion of the sterilization cycle the biological indicator was incubated for seven (7) days and the computer was tested for functionality.

## Data:

**Table 1. Dosimeter Results.**

Cycle	Pass (Y/N)
Anprolene	Y
EOGas	Y

**Table 2. Biological Indicator Results After Seven (7) Days.**

Cycle	Pos/Neg
Anprolene	Neg
EOGas	Neg

**Table 3. Functionality of Computer.**

Cycle	Pass (Y/N)
Anprolene	Y
EOGas	Y

## Conclusion:

Andersen study, PR0804-17: full diagnostics and functionality testing at the conclusion of each sterilization cycle indicated that the laptop's performance had not been impaired.