DOSIMETER® EO Type 5

Name and function of the indicator: indication of the total EO exposure. Conforms to the requirements for Type 5 chemical indicators for EO sterilization processes (ISO 11140-1) (Only for European market).

Device Description
The AN1036 Dosimeter is a single-use chemical indicator for cumulative ethylene oxide exposure. It is an accessory for the EOGas 4® Ethylene Oxide Gas Sterilizer.

The AN1036 Dosimeter contains a proprietary pH indicator in a glass capillary tube that is sealed on one end and mounted on a plastic tray. It is calibrated for a 50°C sterilization temperature, and responds to ethylene oxide concentration and sterilization time. With exposure to ethylene oxide, the indicator turns from yellow-orange to a dark blue color from the open end toward the closed end. The extent of the color change is proportional to cumulative ethylene oxide exposure. The cycle-specific calibration marks represent adequate cumulative ethylene oxide exposure to inactivate a 6-Log Bacillus atrophaeus biological indicator at the location of the AN1036 Dosimeter after 3-hour gas exposures or at the center of endoscope lumens after 6-hour gas exposures. The AN1036 Dosimeter is not a replacement for a biological indicator. The chemical indicator (AN1036 Dosimeter) should be read as soon as practical after the sterilization cycle and the results recorded.

The AN1036 Dosimeter for the EOGas 4 sterilizer has two calibration marks, one for "3-hours EO exposure" and one for "6-hours EO exposure". The depiction below shows a Dosimeter that has not been exposed to EO:

The blue line travels from left to right along the Dosimeter during the cycle. The blue line must reach or pass the cycle-specific calibration mark for the cycle to "pass". In the example below, the line does not reach either calibration mark indicating a "failed" sterilization cycle:

If the blue line reach or passes the "3-hours EO exposure" calibration mark but not the "6-hours EO exposure" calibration mark, a cycle with a 3-hour gas exposure would "pass" whereas a 6-hour gas exposure would "fail":

If the blue line reach or passes the "6-hours EO exposure" calibration mark, both cycle lengths "3-hours EO exposure" or "6-hours EO exposure" would "pass"

Indications for Use
AN1036 Dosimeter is a single-use color change chemical indicator. It is calibrated for sterilization temperature and used to verify adequate cumulative ethylene oxide exposure in a 3-hour or 6-hour gas exposure at 50°C in the EOGas 4 Ethylene Oxide Gas Sterilizer manufactured by Andersen Sterilizers, Inc. Critical process parameters for the cycles are summarized in Table 1.

Table 1. Critical sterilization cycle parameters in the EOGas 4 Ethylene Oxide Gas Sterilizer	
Ethylene Oxide:	17.6 g ± 5%
Temperature:	50°C ± 3°C
Relative Humidity:	35-90%
Ethylene Oxide Exposure Time:	3 hours/6 hours
Total Cycle Time:	3.5 hours/7 hours

Instructions for Use
AN1036 Dosimeters are supplied in sealed foil pouches. Tear at the notch or cut open the foil pouch, remove the Dosimeter, and place one at the center of the load. No safety precautions are required during normal use. After the cycle, confirm that the blue line has reached or passed the cycle-specific calibration mark. This confirms that the time and gas concentration necessary for sterilization have been achieved, and the sterilization cycle is considered to have "passed". If the blue line does not reach the calibration mark, the sterilization cycle is considered to have "failed". Release of sterilized medical devices for use after each sterilization cycle must be based on inactivation of a biological indicator.

Performance
AN1036 Dosimeters comply with the Type 5 performance and labeling requirements of ISO 11140-1. At 54°C, 600 mg/L, 60% RH, and atmospheric pressure, the 3-hour gas exposure calibration mark is reached in 55±5 minutes. At 37°C, 600 mg/L, 60% RH, and atmospheric pressure, the 3-hour gas exposure calibration mark is reached in 165±15 minutes.

Performance Limitations
1) AN1036 Dosimeters are designed for use only in atmospheric pressure EOGas sterilization system;
2) AN1036 Dosimeters are for single use only;
3) The calibration marks are cycle-specific.

Storage Conditions
Store AN1036 Dosimeter sealed in foil pouches at 20-30°C, out of direct sunlight before use. If Dosimeters are stored after use, blue lines are stable (within 15% of the initial length) for 28 days at 20-25°C and for 4 days at 50°C.

Disposal
Dispose of used and unused AN1036 Dosimeters in ordinary trash.

Shelf Life
Expiration date is 3 years from date of manufacture, as indicated on the lot label attached to the foil pouch, if the AN1036 Dosimeter is stored under the recommended conditions.

Manufactured by:
Andersen Sterilizers, Inc.
Health Science Park • 3154 Caroline Dr.
Haw River, NC 27258 USA

Distributed by:
Andersen Products, Inc.
Health Science Park • 3202 Caroline Dr.
Haw River, NC 27258 USA

UK Representative:
H. W. Andersen Products Ltd.
808 Fowler Road
Clacton-On-Sea • Essex CO15 4AA UK

Authorized EU Representative:
H.W. Andersen SARL
12 Place Saint-Hubert
59000 • Lille • France

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