



ANDERSEN
STERILIZERS

ISO 13485:2016

QUALITY MANAGEMENT SYSTEM

Gas Sterilization & Abatement Systems

*Andersen Sterilizers, Inc., Health Science Park,
Haw River, NC 27258 USA*

Sterilizer Accessories Quality Control Certificate

Refill Kit AN#	Lot #	Quantity	Expiration Date
2011.co	257732	108	2030-10-14

Kit Components

☐ AN2005.25	4.74g	EO Gas Filled Ampoule Lot #:		256648	
☐ AN3033.48	7.94g				
☒ AN2023.00	10.5g	EO Gas Filled Ampoule Expiration Date:		2030-10-14	
☐ AN3033.00	17.6g				
☐ N/A					
Lot Avg. Volatile Weight		10.70g	ETO Purity		99.99%
Sterilization Bag	PN#	☐ 822	Lot #	N/A	Expiration Date
		☐ 4646			
		☐ 750004			
		☒ N/A			
Gas Release Bag	PN#	☐ 978	Lot #	N/A	Expiration Date
		☒ N/A			
Dosimeter	AN#	☐ AN87.14	Lot #	N/A	Expiration Date
		☐ AN1087.16			
		☐ AN1036.00			
		☒ N/A			
Humidichips	AN#	☐ AN1071.16	Lot #	N/A	Expiration Date
		☒ N/A			

This product has been manufactured in accordance with the Quality Requirements and Manufacturing Procedures Noted Below.

Product Passed the Following Inspection and Testing:

Packaging DOT	☒
Labeling (EPA/DOT)	☒
In-Process & Final QC Testing	☒

Manufacturing QC:		Date:	2025-12-31
Quality Assurance Review:	Sierra Snow	Date:	2025.12.31